neoBona...

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www.neobona.be
neobona@synlab.be

Références:

- Bianchi DW et al. Genome-wide fetal aneuploidy detection by maternal plasma DNA sequencing. *Obstet Gynecol.* 2012;119:890–901.
- Bianchi DW et al. DNA sequencing versus standard prenatal aneuploidy screening. *N Engl J Med* 2014;370:799–808.
- Cirigliano et al. Performance evaluation and clinical implementation of the NeoBona test, a new paired-end MPSS approach for cfDNA based prenatal screening of common chromosome aneuploidies. *Prenatal Diagnosis* 2016, 36, 23–84.
- Cirigliano et al. First clinical application of paired-end MPSS for cfDNA based prenatal screening of aneuploidies P01.060D ESHG 2016.
- Cirigliano V, Ordonez E, Rueda L, Syngelaki A, Nicolaides KH. Performance evaluation of the NeoBona test, a new paired-end massive parallel shotgun sequencing approach for cfDNA based aneuploidy screening. *Ultrasound in Obstetrics & Gynecology*. DOI:10.1002/uog.17386.
- Futch T et al. Initial clinical laboratory experience in noninvasive prenatal testing for fetal aneuploidy from maternal plasma DNA samples. *Prenat Diagn.* 2013;33:569–574.
- Rava RP et al. Circulating fetal cell-free DNA fractions differ in autosomal aneuploidies and monosomy X. *Clin Chem.* 2014;60:243–250.
- Sehnert AJ et al. Optimal detection of fetal chromosomal abnormalities by massively parallel DNA sequencing of cell-free fetal DNA from maternal blood. *Clin Chem.* 2011;57:1042–1049.
- Srinivasan A et al. Noninvasive detection of fetal subchromosome abnormalities via deep sequencing of maternal plasma. *Am J Human Genet* 2013;92:1–10.
- SYNLAB clinical performance data.

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	SENSIBILITÉ (95% CI)*	SPÉCIFICITÉ (95% CI)*
TRISOMY 21	100% (94.3-100%)	99.96% (99.9-100%)
TRISOMY 18	97.1% (84.7-99.9%)	100% (99.9-100%)
TRISOMY 13	100% (75.3-100%)	99.98% (99.9-100%)

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* Cirigliano V, Ordonez E, Rueda L, Syngelaki A, Nicolaides KH. Performance evaluation of the NeoBona test, a new paired-end massive parallel shotgun sequencing approach for cfDNA based aneuploidy screening. *Ultrasound in Obstetrics & Gynecology*. DOI:10.1002/uog.17386.

* Cirigliano et al. Performance evaluation and clinical implementation of the neoBona test, a new paired-end MPSS approach for cfDNA based prenatal screening of common chromosome aneuploidies. *Prenatal Diagnosis* 2016, 36, 23-84.

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

Cirigliano V, Ordonez E, Rueda L, Syngelaki A, Nicolaides KH. Performance evaluation of the NeoBona test, a new paired-end massive parallel shotgun sequencing approach for cfDNA based aneuploidy screening. *Ultrasound in Obstetrics & Gynecology*. DOI:10.1002/uog.17386.

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